



SPECTRUM PRO – HOSPITAL USE

A BROAD SPECTRUM BACTERICIDE, VIRUCIDE AND FUNGICIDE

Decon Solutions' Spectrum Pro CDD-1 is an EPA registered disinfectant and sanitizer that is effective against a broad spectrum of bacteria, is virucidal, fungicidal, controls odors and inhibits the growth of mold and mildew.

USES - When used as directed Decon Solutions' Spectrum Pro is a proven disinfectant, cleaner, sanitizer, fungicide, mildewstat, and virucide. Spectrum Pro can be applied to walls and floors and other hard surfaces such as tables, chairs, countertops, sinks, glazed tile, and glazed porcelain with a cloth, mop or spray. Can be used on food contact surfaces after a rinse with potable water.

CLEANS AND DECONTAMINATES: viruses such as Hepatitis B Virus, HIV and HBV on precleaned environmental surfaces.

APPROVED HOSPITAL USE: This product meets AOAC Use- Dilution Standards for hospital disinfectants at 400 ppm water hardness.

SANITIZE, DEODORIZE AND CLEAN surfaces quickly and efficiently. Kills common germs and viruses.

Spectrum Pro CDD-1

Hospital Disinfectant & Sanitizer

Bactericidal:

This product is bactericidal according to AOAC Use Dilution Test method on hard inanimate surfaces, modified in the presence of 5% organic serum against:

<i>Brevibacterium ammoniagenes</i>	<i>Proteus vulgaris</i>
<i>Camphylobacter jejuni</i>	<i>Pseudomonas aeruginosa</i>
<i>Corynebacterium ammoniagenes</i>	<i>Pseudomonas aeruginosa (antibiotic resistant)</i>
<i>Enterobacter aerogenes</i>	<i>Salmonella choleraesuis</i>
<i>Enterobacter faecalis</i>	<i>Salmonella marcescens</i>
<i>Enterococcus faecalis</i>	<i>Serratia typhi</i>
<i>Enterococcus faecium (Vancomycin resistant)</i>	<i>Shigella dysenteriae</i>
<i>Escherichia coli</i>	<i>Shigella flexneri</i>
<i>Escherichia coli (Antibiotic resistant)</i>	<i>Shigella sonnei</i>
<i>Escherichia coli 0157:H7</i>	<i>Staphylococcus aureus</i>
<i>Lengionella pneumophila</i>	<i>Staphylococcus aureus (MRSA)</i>
<i>Klebsiella pneumoniae</i>	<i>Staphylococcus epidermidis (antibiotic resistant)</i>
<i>Listeria monocytogenes</i>	<i>Streptococcus phyogenes</i>

Virucidal:

This product was evaluated and found to be effective in the presence of 5% serum and 400 ppm hard water against the following viruses on inanimate hard surface.

<i>Avian influenza A/Turkey/Wisconsin</i>	<i>Infectious Bovine Rhinotracheitis (IBR)</i>
<i>Canine Distemper</i>	<i>Influenza A/Brazil Virus</i>
<i>Herpes Simplex Type 1</i>	<i>Newcastle disease virus</i>
<i>Herpes Simplex Type 2</i>	<i>Porcine Respiratory & Reproductive Virus (PRRSV)</i>
<i>HIV-1 (AIDS virus)</i>	<i>Porcine Rotavirus</i>
<i>Hepatitis C Virus</i>	<i>Pseudorabies virus</i>
<i>Hepatitis B virus</i>	<i>Transmissible Gastroenteritis (TGE)</i>
<i>Human Coronavirus</i>	<i>Vaccinia virus</i>
<i>Canine Coronavirus</i>	<i>Hantavirus</i>

Fungicidal:

Fungicidal against *Candida albicans* and *Trichophyton mentagrophytes* on inanimate surfaces in the presence of 5% blood serum and 400 ppm hard water.

Approved Uses:

- Cross contamination is of major housekeeping concern. This product has been formulated to aid in the reduction of cross-contamination on treated surfaces not only in hospitals, but also in schools, institutions and industry.
- Patient care rooms & facilities, recovery, anesthesia, Emergency Rooms, X-ray cat labs, newborn nurseries, orthopedics, whirlpools, footbaths, respiratory therapy, surgi-centers, labs, blood collection rooms, central supply, housekeeping & janitorial rooms.
- EMS & fire facilities, day care centers and nurseries, sick rooms, acute care institutions, alternate care institutions, and home healthcare institutions.